




Manufacturer & Exporter Of  
Orthopedic Implants



## Orthopedic Implants Catalogue

+923006562639  
+92300 0507475  
orthozenew@yahoo.com  
Shop#1, St#7, Main Bazar Asad Colony,  
Sheikhupura RoadGujranwala  
Pakistan-52250



Manufacture & Exporter of  
Orthopedic Implants



Head Dia : 3mm Hexagonal socket : 1.5mm

Core Dia : 1.2mm Tap Dia : 1.5mm

TPI Dia : 0.40mm Drill Bit : 1.2mm

Shaft Dia : N/A Guide Wire : N/A

Thread Dia : 1.5mm

## 1.5mm Cortical Screw-Self Taping

Length	Stainless Steel 316L	Titanium Grade 5
6mm	ZE-0131S	ZE-0231T
8mm	ZE-0132S	ZE-0232T
10mm	ZE-0133S	ZE-0233T
12mm	ZE-0134S	ZE-0234T
14mm	ZE-0135S	ZE-0235T
16mm	ZE-0136S	ZE-0236T
18mm	ZE-0137S	ZE-0237T
20mm	ZE-0138S	ZE-0238T
22mm	ZE-0139S	ZE-0239T

## 2.0 mm Cortical Screw Self Taping

Head Dia : 4mm Hexagonal socket : 1.5mm


Core Dia : 1.5mm Tap Dia : 2mm

TPI Dia : 0.60mm Drill Bit : 1.5mm

Shaft Dia : N/A Guide Wire : N/A

Thread Dia : 2.0mm

Length	Stainless Steel 316L	Titanium Grade 5
6mm	ZE-0141S	ZE-0241T
8mm	ZE-0142S	ZE-0242T
10mm	ZE-0143S	ZE-0243T
12mm	ZE-0144S	ZE-0244T
14mm	ZE-0145S	ZE-0245T
16mm	ZE-0146S	ZE-0246T
18mm	ZE-0147S	ZE-0247T
20mm	ZE-0148S	ZE-0248T
22mm	ZE-0149S	ZE-0249T



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

## 2.4 mm Cortical Screw Self Taping



Head Dia : 4mm Hexagonal socket : 2.0mm


Core Dia : 1.7mm Tap Dia : 2.4mm

TPI Dia : 0.80mm Drill Bit : 1.8mm

Shaft Dia : N/A Guide Wire : N/A

Thread Dia : 2.4mm

Length	Stainless Steel 316L	Titanium Grade 5
6mm	ZE-0151S	ZE-0251T
8mm	ZE-0152S	ZE-0252T
10mm	ZE-0153S	ZE-0253T
12mm	ZE-0154S	ZE-0254T
14mm	ZE-0155S	ZE-0255T
16mm	ZE-0156S	ZE-0256T
18mm	ZE-0157S	ZE-0257T
20mm	ZE-0158S	ZE-0258T
22mm	ZE-0159S	ZE-0259T
24mm	ZE-0160S	ZE-0260T
26mm	ZE-0161S	ZE-0261T
28mm	ZE-0162S	ZE-0262T
30mm	ZE-0163S	ZE-0263T



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### 2.7 mm Cortical Screw Self Taping



Head Dia : 5mm Hexagonal socket : 2.0mm


Core Dia : 1.9mm Tap Dia : 2.7mm

TPI Dia : 1.0mm Drill Bit : 2.0mm

Shaft Dia : N/A Guide Wire : N/A

Thread Dia: 2.7mm

Length	Stainless Steel 316L	Titanium Grade 5
6mm	ZE-0164S	ZE-0264T
8mm	ZE-0165S	ZE-0265T
10mm	ZE-0166S	ZE-0266T
12mm	ZE-0167S	ZE-0267T
14mm	ZE-0168S	ZE-0268T
16mm	ZE-0169S	ZE-0269T
18mm	ZE-0170S	ZE-0270T
20mm	ZE-0171S	ZE-0271T
22mm	ZE-0172S	ZE-0272T
24mm	ZE-0173S	ZE-0273T
26mm	ZE-0174S	ZE-0274T
28mm	ZE-0175S	ZE-0275T
30mm	ZE-0176S	ZE-0276T
32mm	ZE-0177S	ZE-0277T
34mm	ZE-0178S	ZE-0278T
36mm	ZE-0179S	ZE-0279T
38mm	ZE-0180S	ZE-0280T
40mm	ZE-0181S	ZE-0281T



📍 Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

✉️ orthozenew@yahoo.com




Manufacture & Exporter of  
Orthopedic Implants

3.5 mm Cortical Screw Self Taping



Head Dia : 6mm Hexagonal socket : 2.5mm
Core Dia : 2.4mm Tap Dia : 3.5mm
TPI Dia : 1.25mm Drill Bit : 2.5mm
Shaft Dia : N/A Guide Wire : N/A
Thread Dia: 3.5mm

Length	Stainless Steel 316L	Titanium Grade 5
6mm	ZE-0182S	ZE-0182T
8mm	ZE-0183S	ZE-0183T
10mm	ZE-0184S	ZE-0184T
12mm	ZE-0185S	ZE-0185T
14mm	ZE-0186S	ZE-0186T
16mm	ZE-0187S	ZE-0187T
18mm	ZE-0188S	ZE-0188T
20mm	ZE-0189S	ZE-0189T
22mm	ZE-0190S	ZE-0190T
24mm	ZE-0191S	ZE-0191T
26mm	ZE-0192S	ZE-0192T
28mm	ZE-0193S	ZE-0193T
30mm	ZE-0194S	ZE-0194T
32mm	ZE-0195S	ZE-0195T
34mm	ZE-0196S	ZE-0196T
36mm	ZE-0197S	ZE-0197T
38mm	ZE-0198S	ZE-0198T
40mm	ZE-0199S	ZE-0199T
42mm	ZE-0200S	ZE-0200T
44mm	ZE-0201S	ZE-0201T
46mm	ZE-01202S	ZE-01202T
48mm	ZE-0203S	ZE-0203T
50mm	ZE-0204S	ZE-0204T
52mm	ZE-0205S	ZE-0205T
54mm	ZE-0206S	ZE-0206T
56mm	ZE-0207S	ZE-0207T
58mm	ZE-0208S	ZE-0208T
60mm	ZE-0209S	ZE-0209T
62mm	ZE-0210S	ZE-0210T
64mm	ZE-0211S	ZE-0211T
66mm	ZE-0212S	ZE-0212T
68mm	ZE-0213S	ZE-0213T
70mm	ZE-0214S	ZE-0214T



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com



Manufacture & Exporter of  
Orthopedic Implants

### 4.5 mm Cortical Screw Self Taping



Head Dia : 8mm Hexagonal socket : 3.5mm


Core Dia : 3.0mm Tap Dia : 4.5mm

TPI Dia : 1.75mm Drill Bit : 2.0mm

Shaft Dia : N/A Guide Wire : N/A

Thread Dia: 4.5mm

Length	Stainless Steel 316L	Titanium Grade 5
10mm	ZE-0217S	ZE-0217T
12mm	ZE-0218S	ZE-0218T
14mm	ZE-0219S	ZE-0219T
16mm	ZE-0220S	ZE-0220T
18mm	ZE-0221S	ZE-0221T
20mm	ZE-0222S	ZE-0222T
22mm	ZE-0223S	ZE-0223T
24mm	ZE-0224S	ZE-0224T
26mm	ZE-0225S	ZE-0225T
28mm	ZE-0226S	ZE-0226T
30mm	ZE-0227S	ZE-0227T
32mm	ZE-0228S	ZE-0228T
34mm	ZE-0229S	ZE-0229
36mm	ZE-0230S	ZE-0230T
38mm	ZE-0231S	ZE-0231T
40mm	ZE-0232S	ZE-0232T
42mm	ZE-0233S	ZE-0233T
44mm	ZE-0234S	ZE-0234T
46mm	ZE-0235S	ZE-0235T
48mm	ZE-0236S	ZE-0236T
50mm	ZE-0237S	ZE-0237T
52mm	ZE-0238S	ZE-0238T
54mm	ZE-0239S	ZE-0239T
56mm	ZE-0240S	ZE-0240T
58mm	ZE-0241S	ZE-0241T
60mm	ZE-0242S	ZE-0242T
62mm	ZE-0243S	ZE-0243T
64mm	ZE-0244S	ZE-0244T
66mm	ZE-0245S	ZE-0245T
68mm	ZE-0246S	ZE-0246T
70mm	ZE-0247S	ZE-0247T




Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

1.5 mm Safety Locking Screws (LCP)

Length	Stainless Steel 316L	Titanium Grade 5
6mm	ZE-0248S	ZE-0248T
8mm	ZE-0249S	ZE-0249T
10mm	ZE-0250S	ZE-0250T
12mm	ZE-0251S	ZE-0251T
14mm	ZE-0252S	ZE-0252T
16mm	ZE-0253S	ZE-0253T
18mm	ZE-0254S	ZE-0254T
20mm	ZE-0255S	ZE-0255T
22mm	ZE-0256S	ZE-0256T
24mm	ZE-0257S	ZE-0257T
26mm	ZE-0258S	ZE-0258T
28mm	ZE-0259S	ZE-0259T
30mm	ZE-0260S	ZE-0260T



📍 Shop No1, St#7, Main Bazar Asad Colony, Sheikhupura Road  
Gujranwala Pakistan-52250.

✉️ orthozenew@yahoo.com




Manufacture & Exporter of  
Orthopedic Implants

### 2.0 mm Safety Locking Screw (LCP)



Length	Stainless Steel 316L	Titanium Grade 5
6mm	ZE-0261S	ZE-0261T
8mm	ZE-0262S	ZE-0262T
10mm	ZE-0263S	ZE-0263T
12mm	ZE-0264S	ZE-0264T
14mm	ZE-0265S	ZE-0265T
16mm	ZE-0266S	ZE-0266T
18mm	ZE-0267S	ZE-0267T
20mm	ZE-0268S	ZE-0268T
22mm	ZE-0269S	ZE-0269T
24mm	ZE-0270S	ZE-0270T
26mm	ZE-0271S	ZE-0271T
28mm	ZE-0272S	ZE-0272T
30mm	ZE-0273S	ZE-0273T



📍 Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

✉️ [orthozenew@yahoo.com](mailto:orthozenew@yahoo.com)



### 2.4 mm Safety Locking Screws (LCP)



Head Dia : 3.4mm Hexagonal socket : 2mm


Core Dia : 1.9mm Tap Dia : N/A

TPI Dia : 0.6mm Drill Bit : 1.8mm

Shaft Dia : N/A Guide Wire : N/A

Thread Dia: 2.4mm

Length	Stainless Steel 316L	Titanium Grade 5
8mm	ZE-0274S	ZE-0274T
10mm	ZE-0275S	ZE-0275T
12mm	ZE-0276S	ZE-0276T
14mm	ZE-0277S	ZE-0277T
16mm	ZE-0278S	ZE-0278T
18mm	ZE-0279S	ZE-0279T
20mm	ZE-0280S	ZE-0280T
22mm	ZE-0281S	ZE-0281T
24mm	ZE-0282S	ZE-0282T
26mm	ZE-0283S	ZE-0283T
28mm	ZE-0284S	ZE-0284T
30mm	ZE-0285S	ZE-0285T
32mm	ZE-0286S	ZE-0286T
34mm	ZE-0287S	ZE-0287T
36mm	ZE-0288S	ZE-0288T
38mm	ZE-0289S	ZE-0289T
40mm	ZE-0290S	ZE-0290T



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.


orthozenew@yahoo.com

### 2.7 mm Safety Locking Screw (LCP)



Length	Stainless Steel 316L	Titanium Grade 5
8mm	ZE-0291S	ZE-0291T
10mm	ZE-0292S	ZE-0292T
12mm	ZE-0293S	ZE-0293T
14mm	ZE-0294S	ZE-0294T
16mm	ZE-0295S	ZE-0295T
18mm	ZE-0296S	ZE-0296T
20mm	ZE-0297S	ZE-0297T
22mm	ZE-0298S	ZE-0298T
24mm	ZE-0299S	ZE-0299T
26mm	ZE-0300S	ZE-0300T
28mm	ZE-0301S	ZE-0301T
30mm	ZE-0302S	ZE-0302T
32mm	ZE-0303S	ZE-0303T
34mm	ZE-0304S	ZE-0304T
36mm	ZE-0305S	ZE-0305T
38mm	ZE-0306S	ZE-0306T
40mm	ZE-0307S	ZE-0307T

Head Dia : 3.4mm Hexagonal socket : 2mm
Core Dia : 2.1mm Tap Dia : N/A
TPI Dia : 0.6mm Drill Bit : 2.0mm
Shaft Dia : N/A Guide Wire : N/A
Thread Dia: 2.7mm



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### 3.5 mm Safety Locking Screws (LCP)



Head Dia : 3.4mm Hexagonal socket : 2mm


Core Dia : 1.9mm Tap Dia : N/A

TPI Dia : 0.6mm Drill Bit : 1.8mm

Shaft Dia : N/A Guide Wire : N/A

Thread Dia: 2.4mm


Length	Stainless Steel 316L	Titanium Grade 5
10mm	ZE-0308S	ZE-0308T
12mm	ZE-0309S	ZE-0309T
14mm	ZE-0310S	ZE-0310T
16mm	ZE-0311S	ZE-0311T
18mm	ZE-0312S	ZE-0312T
20mm	ZE-0313S	ZE-0313T
22mm	ZE-0314S	ZE-0314T
24mm	ZE-0315S	ZE-0315T
26mm	ZE-0316S	ZE-0316T
28mm	ZE-0317S	ZE-0317T
30mm	ZE-0318S	ZE-0318T
32mm	ZE-0319S	ZE-0319T
34mm	ZE-0320S	ZE-0320T
36mm	ZE-0321S	ZE-0321T
38mm	ZE-0322S	ZE-0322T
40mm	ZE-0323S	ZE-0323T
42mm	ZE-0324S	ZE-0324T
44mm	ZE-0325S	ZE-0325T
46mm	ZE-0326S	ZE-0326T
48mm	ZE-0327S	ZE-0327T
50mm	ZE-0328S	ZE-0328T
52mm	ZE-0329S	ZE-0329T
54mm	ZE-0330S	ZE-0330T
56mm	ZE-0331S	ZE-0331T
60mm	ZE-0332S	ZE-0332T
65mm	ZE-0333S	ZE-0333T
70mm	ZE-0334S	ZE-0334T
75mm	ZE-0335S	ZE-0335T
80mm	ZE-0336S	ZE-0336T
85mm	ZE-0337S	ZE-0337T
90mm	ZE-0338S	ZE-0338T



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### 5 mm Safety Locking Screw (LCP)



Length	Stainless Steel 316L	Titanium Grade 5
10mm	ZE-0339S	ZE-0339T
12mm	ZE-0340S	ZE-0340T
14mm	ZE-0341S	ZE-0341T
16mm	ZE-0342S	ZE-0342T
18mm	ZE-0343S	ZE-0343T
20mm	ZE-0344S	ZE-0344T
22mm	ZE-0345S	ZE-0345T
24mm	ZE-0346S	ZE-0346T
26mm	ZE-0347S	ZE-0347T
28mm	ZE-0348S	ZE-0348T
30mm	ZE-0349S	ZE-0349T
32mm	ZE-0350S	ZE-0350T
34mm	ZE-0351S	ZE-0351T
36mm	ZE-0352S	ZE-0352T
38mm	ZE-0353S	ZE-0353T
40mm	ZE-0354S	ZE-0354T
42mm	ZE-0355S	ZE-0355T
44mm	ZE-0356S	ZE-0356T
46mm	ZE-0357S	ZE-0357T
48mm	ZE-0358S	ZE-0358T
50mm	ZE-0359S	ZE-0359T
52mm	ZE-0360S	ZE-0360T
54mm	ZE-0361S	ZE-0361T
56mm	ZE-0362S	ZE-0362T
60mm	ZE-0363S	ZE-0363T
65mm	ZE-0364S	ZE-0364T
70mm	ZE-0365S	ZE-0365T
75mm	ZE-0366S	ZE-0366T
80mm	ZE-0367S	ZE-0367T
85mm	ZE-0368S	ZE-0368T
90mm	ZE-0369S	ZE-0369T

Head Dia : 6.5mm Hexagonal socket : 3.5mm


Core Dia : 4.44mm Tap Dia : N/A

TPI Dia : 1.0mm Drill Bit : 4.3mm

Shaft Dia : N/A Guide Wire : N/A

Thread Dia: 5mm

Manufacture & Exporter of  
Orthopedic Implants



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.


orthozenew@yahoo.com

### 3.5 mm Cancellous Screw Partial Thread



Head Dia : 6 mm Hexagonal socket : 2 .5mm
Core Dia : 1.9 mm Tap Dia : N/A
TPI Dia : 1.75mm Drill Bit : 2 mm
Shaft Dia : N/A Guide Wire : N/A
Thread Dia: 3.5mm

Length	Stainless Steel 316L	Titanium Grade 5
10mm	ZE-0370S	ZE-0370T
12mm	ZE-0371S	ZE-0371T
14mm	ZE-0372S	ZE-0372T
16mm	ZE-0373S	ZE-0373T
18mm	ZE-0374S	ZE-0374T
20mm	ZE-0375S	ZE-0375T
22mm	ZE-0376S	ZE-0376T
24mm	ZE-0377S	ZE-0377T
26mm	ZE-0378S	ZE-0378T
28mm	ZE-0379S	ZE-0379T
30mm	ZE-0380S	ZE-0380T
35mm	ZE-0381S	ZE-0381T
40mm	ZE-0382S	ZE-0382T
45mm	ZE-0383S	ZE-0383T
50mm	ZE-0384S	ZE-0384T
55mm	ZE-0385S	ZE-0385T
60mm	ZE-0386S	ZE-0386T
65mm	ZE-0387S	ZE-0387T
70mm	ZE-0388S	ZE-0388T




Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### 3.5 mm Cancellous Screw Full Thread

Length	Stainless Steel 316L	Titanium Grade 5
10mm	ZE-0389S	ZE-0389T
12mm	ZE-0390S	ZE-0390T
14mm	ZE-0391S	ZE-0391T
16mm	ZE-0392S	ZE-0392T
18mm	ZE-0393S	ZE-0393T
20mm	ZE-0394S	ZE-0394T
22mm	ZE-0395S	ZE-0395T
24mm	ZE-0396S	ZE-0396T
26mm	ZE-0397S	ZE-0397T
28mm	ZE-0398S	ZE-0398T
30mm	ZE-0399S	ZE-0399T
35mm	ZE-0400S	ZE-0400T
40mm	ZE-0401S	ZE-0401T
45mm	ZE-0402S	ZE-0402T
50mm	ZE-0403S	ZE-0403T
55mm	ZE-0404S	ZE-0404T
60mm	ZE-0405S	ZE-0405T
65mm	ZE-0406S	ZE-0406T
70mm	ZE-0407S	ZE-0407T



Head Dia : 6 mm Hexagonal socket : 2 .5mm


Core Dia : 1.9 mm Tap Dia : 3.5mm

TPI Dia : 1.75mm Drill Bit : 2 mm

Shaft Dia : N/A Guide Wire : N/A

Thread Dia: 3.5mm

Manufacture & Exporter of  
Orthopedic Implants



📍 Shop No1, St#7, Main Bazar Asad Colony, Sheikhupura Road  
Gujranwala Pakistan-52250.

✉️ orthozenew@yahoo.com

### 4.0 mm Cancellous Screw Partial Thread



Head Dia : 6 mm Hexagonal socket : 2 .5mm


Core Dia : 1.9 mm Tap Dia : N/A

TPI Dia : 1.75mm Drill Bit : 2 mm

Shaft Dia : 2.5mm Guide Wire : N/A

Thread Dia: 4mm

Length	Stainless Steel 316L	Titanium Grade 5
10mm	ZE-0408S	ZE-0408T
12mm	ZE-0409S	ZE-0409T
14mm	ZE-0410S	ZE-0410T
16mm	ZE-0411S	ZE-0411T
18mm	ZE-0412S	ZE-0412T
20mm	ZE-0413S	ZE-0413T
22mm	ZE-0414S	ZE-0414T
24mm	ZE-0415S	ZE-0415T
26mm	ZE-0416S	ZE-0416T
28mm	ZE-0417S	ZE-0417T
30mm	ZE-0418S	ZE-0418T
35mm	ZE-0419S	ZE-0419T
40mm	ZE-0420S	ZE-0420T
45mm	ZE-0421S	ZE-0421T
50mm	ZE-0422S	ZE-0422T
55mm	ZE-0423S	ZE-0423T
60mm	ZE-0424S	ZE-0424T
65mm	ZE-0425S	ZE-0425T
70mm	ZE-0426S	ZE-0426T



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### 4 mm Cancellous Screw Full Thread

Length	Stainless Steel 316L	Titanium Grade 5
10mm	ZE-0427S	ZE-0427T
12mm	ZE-0428S	ZE-0428T
14mm	ZE-0429S	ZE-0429T
16mm	ZE-0430S	ZE-0430T
18mm	ZE-0431S	ZE-0431T
20mm	ZE-0432S	ZE-0432T
22mm	ZE-0433S	ZE-0433T
24mm	ZE-0434S	ZE-0434T
26mm	ZE-0435S	ZE-0435T
28mm	ZE-0436S	ZE-0436T
30mm	ZE-0437S	ZE-0437T
35mm	ZE-0438S	ZE-0438T
40mm	ZE-0439S	ZE-0439T
45mm	ZE-0440S	ZE-0440T
50mm	ZE-0441S	ZE-0441T
55mm	ZE-0442S	ZE-0442T
60mm	ZE-0443S	ZE-0443T
65mm	ZE-0444S	ZE-0444T
70mm	ZE-0445S	ZE-0445T




Head Dia : 6 mm Hexagonal socket : 2 .5mm

Core Dia : 1.9 mm Tap Dia : 4 mm

TPI Dia : 1.75mm Drill Bit : 2 mm

Shaft Dia : N/A Guide Wire : N/A

Thread Dia: 4mm



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com



### 6.5 mm Cancellous Screw-16mm Thread

Length	Stainless Steel 316L	Titanium Grade 5
25mm	ZE-0446S	ZE-0446T
30mm	ZE-0447S	ZE-0447T
35mm	ZE-0448S	ZE-0448T
40mm	ZE-0449S	ZE-0449T
45mm	ZE-0450S	ZE-0450T
50mm	ZE-0451S	ZE-0451T
55mm	ZE-0452S	ZE-0452T
60mm	ZE-0453S	ZE-0453T
65mm	ZE-0454S	ZE-0454T
70mm	ZE-0455S	ZE-0455T
75mm	ZE-0456S	ZE-0456T
80mm	ZE-0457S	ZE-0457T
85mm	ZE-0458S	ZE-0458T
90mm	ZE-0459S	ZE-0459T
95mm	ZE-0460S	ZE-0460T
100mm	ZE-0461S	ZE-0461T
105mm	ZE-0462S	ZE-0462T
110mm	ZE-0463S	ZE-0463T




Head Dia : 8 mm Hexagonal socket : 3.5mm

Core Dia : 3 mm Tap Dia : N/A

TPI Dia : 2.8mm Drill Bit : 3.2 mm

Shaft Dia : 4.5 Guide Wire : N/A

Thread Dia: 6.5mm



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### 6.5 mm Cancellous Screw-32mm Thread

Head Dia : 8 mm Hexagonal socket : 3.5mm

Core Dia : 3 mm Tap Dia : N/A


TPI Dia : 2.8mm Drill Bit : 3.2 mm

Shaft Dia : 4.5 Guide Wire : N/A

Thread Dia: 6.5mm




Length	Stainless Steel 316L	Titanium Grade 5
30mm	ZE-0464S	ZE-0464T
35mm	ZE-0465S	ZE-0465T
40mm	ZE-0466S	ZE-0466T
45mm	ZE-0467S	ZE-0467T
50mm	ZE-0468S	ZE-0468T
55mm	ZE-0469S	ZE-0469T
60mm	ZE-0470S	ZE-0470T
65mm	ZE-0471S	ZE-0471T
70mm	ZE-0472S	ZE-0472T
75mm	ZE-0473S	ZE-0473T
80mm	ZE-0474S	ZE-0474T
85mm	ZE-0475S	ZE-0475T
90mm	ZE-0476S	ZE-0476T
95mm	ZE-0477S	ZE-0477T
100mm	ZE-0478S	ZE-0478T
105mm	ZE-0479S	ZE-0479T
110mm	ZE-0480S	ZE-0480T



📍 Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

✉️ orthozenew@yahoo.com

### 6.5 mm Cancellous Screw-Full Thread



Head Dia : 8 mm Hexagonal socket : 3.5mm


Core Dia : 3 mm Tap Dia : N/A

TPI Dia : 2.8mm Drill Bit : 3.2 mm

Shaft Dia : 4.5 Guide Wire : N/A

Thread Dia: 6.5mm

Length	Stainless Steel 316L	Titanium Grade 5
30mm	ZE-0481S	ZE-0481T
35mm	ZE-0482S	ZE-0482T
40mm	ZE-0483S	ZE-0483T
45mm	ZE-0484S	ZE-0484T
50mm	ZE-0485S	ZE-0485T
55mm	ZE-0486S	ZE-0486T
60mm	ZE-0487S	ZE-0487T
65mm	ZE-0488S	ZE-0488T
70mm	ZE-0489S	ZE-0489T
75mm	ZE-0490S	ZE-0490T
80mm	ZE-0491S	ZE-0491T
85mm	ZE-0492S	ZE-0492T
90mm	ZE-0493S	ZE-0493T
95mm	ZE-0494S	ZE-0494T
100mm	ZE-0495S	ZE-0495T
105mm	ZE-0496S	ZE-0496T
110mm	ZE-0497S	ZE-0497T



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### 4 mm LCP Cancellous Screw-Partial Thread

Head Dia : 5 mm Hexagonal socket : 2.5mm

Core Dia : 2.4 mm Tap Dia : N/A

TPI Dia : 1.75mm Drill Bit : 2.5 mm

Shaft Dia : 2.5 Guide Wire : N/A

Thread Dia: 4mm



Length	Stainless Steel 316L	Titanium Grade 5
10mm	ZE-0498S	ZE-0498T
12mm	ZE-0499S	ZE-0499T
14mm	ZE-0500S	ZE-0500T
16mm	ZE-0501S	ZE-0501T
18mm	ZE-0502S	ZE-0502T
20mm	ZE-0503S	ZE-0503T
22mm	ZE-0504S	ZE-0504T
24mm	ZE-0505S	ZE-0505T
26mm	ZE-0506S	ZE-0506T
28mm	ZE-0507S	ZE-0507T
30mm	ZE-0508S	ZE-0508T
35mm	ZE-0509S	ZE-0509T
40mm	ZE-0510S	ZE-0510T
45mm	ZE-0511S	ZE-0511T
50mm	ZE-0512S	ZE-0512T
55mm	ZE-0513S	ZE-0513T
60mm	ZE-0514S	ZE-0514T
65mm	ZE-0515S	ZE-0515T
70mm	ZE-0516S	ZE-0516T


Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com



Manufacture & Exporter of  
Orthopedic Implants

### 4 mm LCP Cancellous Screw-Full Thread



Head Dia : 5 mm Hexagonal socket : 2.5mm


Core Dia : 2.4 mm Tap Dia : N/A

TPI Dia : 1.75mm Drill Bit : 2.5 mm

Shaft Dia : N/A Guide Wire : N/A

Thread Dia: 4mm

Length	Stainless Steel 316L	Titanium Grade 5
10mm	ZE-0517S	ZE-0517T
12mm	ZE-0518S	ZE-0518T
14mm	ZE-0519S	ZE-0519T
16mm	ZE-0520S	ZE-0520T
18mm	ZE-0521S	ZE-0521T
20mm	ZE-0522S	ZE-0522T
22mm	ZE-0523S	ZE-0523T
24mm	ZE-0524S	ZE-0524T
26mm	ZE-0525S	ZE-0525T
28mm	ZE-0526S	ZE-0526T
30mm	ZE-0527S	ZE-0527T
35mm	ZE-0528S	ZE-0528T
40mm	ZE-0529S	ZE-0529T
45mm	ZE-0530S	ZE-0530T
50mm	ZE-0531S	ZE-0531T
55mm	ZE-0532S	ZE-0532T
60mm	ZE-0533S	ZE-0533T
65mm	ZE-0534S	ZE-0534T
70mm	ZE-0535S	ZE-0535T



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.


orthozenew@yahoo.com

**5 mm LCP Cancellous Screw-32mm Thread**



Head Dia :6.5 mm Hexagonal socket : 3.5mm
Core Dia : 3.4 mm Tap Dia : N/A
TPI Dia : 2.8mm Drill Bit : 3.5 mm
Shaft Dia : 4.5mm Guide Wire : N/A
Thread Dia: 5mm

Length	Stainless Steel 316L	Titanium Grade 5
30mm	ZE-0536S	ZE-0536T
35mm	ZE-0537S	ZE-0537T
40mm	ZE-0538S	ZE-0538T
45mm	ZE-0539S	ZE-0539T
50mm	ZE-0540S	ZE-0540T
55mm	ZE-0541S	ZE-0541T
60mm	ZE-0542S	ZE-0542T
65mm	ZE-0543S	ZE-0543T
70mm	ZE-0544S	ZE-0544T
75mm	ZE-0545S	ZE-0545T
80mm	ZE-0546S	ZE-0546T
85mm	ZE-0547S	ZE-0547T
90mm	ZE-0548S	ZE-0548T
95mm	ZE-0549S	ZE-0549T
100mm	ZE-0550S	ZE-0550T
105mm	ZE-0551S	ZE-0551T
110mm	ZE-0552S	ZE-0552T




Shop No 1, Sit 7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### 5 mm LCP Cancellous Screw-Full Thread

Length	Stainless Steel 316L	Titanium Grade 5
18mm	ZE-0553S	ZE-0553T
20mm	ZE-0554S	ZE-0554T
22mm	ZE-0555S	ZE-0555T
24mm	ZE-0556S	ZE-0556T
26mm	ZE-0557S	ZE-0557T
28mm	ZE-0558S	ZE-0558T
30mm	ZE-0559S	ZE-0559T
32mm	ZE-0560S	ZE-0560T
34mm	ZE-0561S	ZE-0561T
36mm	ZE-0562S	ZE-0562T
38mm	ZE-0563S	ZE-0563T
40mm	ZE-0564S	ZE-0564T
45mm	ZE-0565S	ZE-0565T
50mm	ZE-0566S	ZE-0566T
55mm	ZE-0567S	ZE-0567T
60mm	ZE-0568S	ZE-0568T
65mm	ZE-0569S	ZE-0569T
70mm	ZE-0570S	ZE-0570T
75mm	ZE-0571S	ZE-0571T
80mm	ZE-0572S	ZE-0572T
85mm	ZE-0573S	ZE-0573T
90mm	ZE-0574S	ZE-0574T
95mm	ZE-0575S	ZE-0575T
100mm	ZE-0576S	ZE-0576T
105mm	ZE-0577S	ZE-0577T
110mm	ZE-0578S	ZE-0578T




Head Dia :6.5 mm Hexagonal socket : 3.5mm

Core Dia : 3.4 mm Tap Dia : N/A

TPI Dia : 2.8mm Drill Bit : 3.5 mm

Shaft Dia : 4.5mm Guide Wire : N/A

Thread Dia: 5mm



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### 6.5 mm LCP Cancellous Screw-16mm Thread

Length	Stainless Steel 316L	Titanium Grade 5
20mm	ZE-0579S	ZE-0579T
22mm	ZE-0580S	ZE-0580T
24mm	ZE-0581S	ZE-0581T
26mm	ZE-0582S	ZE-0582T
28mm	ZE-0583S	ZE-0583T
30mm	ZE-0584S	ZE-0584T
32mm	ZE-0585S	ZE-0585T
34mm	ZE-0586S	ZE-0586T
36mm	ZE-0587S	ZE-0587T
38mm	ZE-0588S	ZE-0588T
40mm	ZE-0589S	ZE-0589T
45mm	ZE-0590S	ZE-0590T
50mm	ZE-0591S	ZE-0591T
55mm	ZE-0592S	ZE-0592T
60mm	ZE-0593S	ZE-0593T
65mm	ZE-0584S	ZE-0584T
70mm	ZE-0595S	ZE-0595T
75mm	ZE-0596S	ZE-0596T
80mm	ZE-0597S	ZE-0597T
85mm	ZE-0598S	ZE-0598T
90mm	ZE-0599S	ZE-0599T
95mm	ZE-0600S	ZE-0600T
100mm	ZE-0602S	ZE-0602T
105mm	ZE-0602S	ZE-0602T
110mm	ZE-0603S	ZE-0603T




Head Dia :8 mm Hexagonal socket : 3.5mm

Core Dia : 3.4 mm Tap Dia : N/A

TPI Dia : 2.8mm Drill Bit : 3.5 mm

Shaft Dia : 4.5mm Guide Wire : N/A

Thread Dia: 6.5mm



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com



### 6.5 mm LCP Cancellous Screw-32mm Thread

Length	Stainless Steel 316L	Titanium Grade 5
30mm	ZE-0604S	ZE-0604T
32mm	ZE-0605S	ZE-0605T
34mm	ZE-0606S	ZE-0606T
36mm	ZE-0607S	ZE-0607T
38mm	ZE-0608S	ZE-0608T
40mm	ZE-0609S	ZE-0609T
45mm	ZE-0610S	ZE-0610T
50mm	ZE-0611S	ZE-0611T
55mm	ZE-0612S	ZE-0612T
60mm	ZE-0613S	ZE-0613T
65mm	ZE-0614S	ZE-0614T
70mm	ZE-0615S	ZE-0615T
75mm	ZE-0616S	ZE-0616T
80mm	ZE-0617S	ZE-0617T
85mm	ZE-0618S	ZE-0618T
90mm	ZE-0619S	ZE-0619T
95mm	ZE-0620S	ZE-0620T
100mm	ZE-0621S	ZE-0621T




Head Dia :8 mm Hexagonal socket : 3.5mm

Core Dia : 3.4 mm Tap Dia : N/A

TPI Dia : 2.8mm Drill Bit : 3.5 mm

Shaft Dia : 4.5mm Guide Wire : N/A

Thread Dia: 6.5mm



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Roa  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### 6.5 mm LCP Cancellous Screw Full Thread

Length	Stainless Steel 316L	Titanium Grade 5
20mm	ZE-0622S	ZE-0622T
22mm	ZE-0623S	ZE-0623T
24mm	ZE-0624S	ZE-0624T
26mm	ZE-0625S	ZE-0625T
28mm	ZE-0626S	ZE-0626T
30mm	ZE-0627S	ZE-0627T
32mm	ZE-0628S	ZE-0628T
34mm	ZE-0629S	ZE-0629T
36mm	ZE-0630S	ZE-0630T
38mm	ZE-0631S	ZE-0631T
40mm	ZE-0632S	ZE-0632T
45mm	ZE-0633S	ZE-0633T
50mm	ZE-0634S	ZE-0634T
55mm	ZE-0635S	ZE-0635T
60mm	ZE-0636S	ZE-0636T
65mm	ZE-0637S	ZE-0637T
70mm	ZE-0638S	ZE-0638T
75mm	ZE-0639S	ZE-0639T
80mm	ZE-0640S	ZE-0640T
85mm	ZE-0641S	ZE-0641T
90mm	ZE-0642S	ZE-0642T
95mm	ZE-0643S	ZE-0643T
100mm	ZE-0644S	ZE-0644T
105mm	ZE-0645S	ZE-0645T
110mm	ZE-0646S	ZE-0646T




Head Dia : 8 mm Hexagonal socket : 3.5mm

Core Dia : 3.4 mm Tap Dia : N/A

TPI Dia : 2.8mm Drill Bit : 3.5 mm

Shaft Dia : 4.5mm Guide Wire : N/A

Thread Dia : 6.5mm



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### 4 mm Cannulated Cancellous Screw Partial Thread

Length	Stainless Steel 316L	Titanium Grade 5
20mm	ZE-0647S	ZE-0647T
22mm	ZE-0648S	ZE-0648T
24mm	ZE-0649S	ZE-0649T
26mm	ZE-0650S	ZE-0650T
28mm	ZE-0651S	ZE-0651T
30mm	ZE-0652S	ZE-0652T
32mm	ZE-0653S	ZE-0653T
34mm	ZE-0654S	ZE-0654T
36mm	ZE-0655S	ZE-0655T
38mm	ZE-0656S	ZE-0656T
40mm	ZE-0657S	ZE-0657T
45mm	ZE-0658S	ZE-0658T
50mm	ZE-0659S	ZE-0659T
55mm	ZE-0660S	ZE-0660T
60mm	ZE-0661S	ZE-0661T
65mm	ZE-0662S	ZE-0662T
70mm	ZE-0663S	ZE-0663T
75mm	ZE-0664S	ZE-0664T
80mm	ZE-0665S	ZE-0665T




Head Dia : 6 mm Hexagonal socket : 3.5mm

Core Dia : 2.6 mm Tap Dia : N/A

TPI Dia : 1.75 mm Drill Bit : 2.7 mm

Shaft Dia : 2.7mm Guide Wire : 1.2mm

Thread Dia: 4 mm



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### 4 mm Cannulated Cancellous Screw Full Thread

Length	Stainless Steel 316L	Titanium Grade 5
20mm	ZE-0666S	ZE-0666T
22mm	ZE-0667S	ZE-0667T
24mm	ZE-0668S	ZE-0668T
26mm	ZE-0669S	ZE-0669T
28mm	ZE-0670S	ZE-0670T
30mm	ZE-0671S	ZE-0671T
32mm	ZE-0672S	ZE-0672T
34mm	ZE-0673S	ZE-0673T
36mm	ZE-0674S	ZE-0674T
38mm	ZE-0675S	ZE-0675T
40mm	ZE-0676S	ZE-0676T
45mm	ZE-0677S	ZE-0677T
50mm	ZE-0678S	ZE-0678T
55mm	ZE-0679S	ZE-0679T
60mm	ZE-0680S	ZE-0680T
65mm	ZE-0681S	ZE-0681T
70mm	ZE-0682S	ZE-0682T
75mm	ZE-0683S	ZE-0683T
80mm	ZE-0684S	ZE-0684T




Head Dia : 6 mm Hexagonal socket : 3.5mm

Core Dia : 2.6 mm Tap Dia : N/A

TPI Dia : 1.75 mm Drill Bit : 2.7 mm

Shaft Dia : 2.7mm Guide Wire : 1.2mm

Thread Dia: 4 mm



📍 Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

✉️ orthozenew@yahoo.com

### 6.5 mm Cannulated Cancellous Screw 16mm Thread

Length	Stainless Steel 316L	Titanium Grade 5
25mm	ZE-0685S	ZE-0685T
30mm	ZE-0686S	ZE-0686T
35mm	ZE-0687S	ZE-0687T
40mm	ZE-0688S	ZE-0688T
45mm	ZE-0689S	ZE-0689T
50mm	ZE-0690S	ZE-0690T
55mm	ZE-0691S	ZE-0691T
60mm	ZE-0692S	ZE-0692T
65mm	ZE-0693S	ZE-0693T
70mm	ZE-0694S	ZE-0694T
75mm	ZE-0695S	ZE-0695T
80mm	ZE-0696S	ZE-0696T
85mm	ZE-0697S	ZE-0697T
90mm	ZE-0698S	ZE-0698T
95mm	ZE-0699S	ZE-0699T
100mm	ZE-0700S	ZE-0700T
105mm	ZE-0701S	ZE-0701T
110mm	ZE-0702S	ZE-0702T




Head Dia : 8 mm Hexagonal socket : 3.5mm

Core Dia : 4.4 mm Tap Dia : N/A

TPI Dia : 2.8 mm Drill Bit : 4.5 mm

Shaft Dia : 4.5mm Guide Wire : 1.8mm

Thread Dia: 6.5 mm



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### 6.5 mm Cannulated Cancellous Screw 32mm Thread

Length	Stainless Steel 316L	Titanium Grade 5
30mm	ZE-0703S	ZE-0703T
35mm	ZE-0704S	ZE-0704T
40mm	ZE-0705S	ZE-0705T
45mm	ZE-0706S	ZE-0706T
50mm	ZE-0707S	ZE-0707T
55mm	ZE-0708S	ZE-0708T
60mm	ZE-0709S	ZE-0709T
65mm	ZE-0710S	ZE-0710T
70mm	ZE-0711S	ZE-0711T
75mm	ZE-0712S	ZE-0712T
80mm	ZE-0713S	ZE-0713T
85mm	ZE-0714S	ZE-0714T
90mm	ZE-0715S	ZE-0715T
95mm	ZE-0716S	ZE-0716T
100mm	ZE-0717S	ZE-0717T
105mm	ZE-0718S	ZE-0718T
110mm	ZE-0719S	ZE-0719T




Head Dia :8 mm Hexagonal socket : 3.5mm

Core Dia : 4.4 mm Tap Dia : N/A

TPI Dia : 2.8 mm Drill Bit : 4.5 mm

Shaft Dia : 4.5mm Guide Wire : 1.8mm

Thread Dia: 6.5 mm



Shop No1, St#7, Main Bazar Asad Colony, Sheikhupura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### 6.5 mm Cannulated Cancellous Screw Full Thread

Length	Stainless Steel 316L	Titanium Grade 5
30mm	ZE-0720S	ZE-0720T
35mm	ZE-0721S	ZE-0721T
40mm	ZE-0722S	ZE-0722T
45mm	ZE-0723S	ZE-0723T
50mm	ZE-0724S	ZE-0724T
55mm	ZE-0725S	ZE-0725T
60mm	ZE-0726S	ZE-0726T
65mm	ZE-0727S	ZE-0727T
70mm	ZE-0728S	ZE-0728T
75mm	ZE-0729S	ZE-0729T
80mm	ZE-0730S	ZE-0730T
85mm	ZE-0731S	ZE-0731T
90mm	ZE-0732S	ZE-0732T
95mm	ZE-0733S	ZE-0733T
100mm	ZE-0734S	ZE-0734T
105mm	ZE-0735S	ZE-0735T
110mm	ZE-0736S	ZE-0736T




Head Dia :8 mm Hexagonal socket : 3.5mm

Core Dia : 4.4 mm Tap Dia : N/A

TPI Dia : 2.8 mm Drill Bit : 4.5 mm

Shaft Dia : 4.5mm Guide Wire : 1.8mm

Thread Dia: 6.5 mm



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### 4.5 mm Malleolar Screw Self Tapping

Length	Stainless Steel 316L	Titanium Grade 5
25mm	ZE-0737S	ZE-0737T
30mm	ZE-0738S	ZE-0738T
35mm	ZE-0739S	ZE-0739T
40mm	ZE-0740S	ZE-0740T
45mm	ZE-0741S	ZE-0741T
50mm	ZE-0742S	ZE-0742T
55mm	ZE-0743S	ZE-0743T
60mm	ZE-0744S	ZE-0744T
65mm	ZE-0745S	ZE-0745T
70mm	ZE-0746S	ZE-0746T
75mm	ZE-0747S	ZE-0747T
80mm	ZE-0748S	ZE-0748T
85mm	ZE-0749S	ZE-0749T
90mm	ZE-0750S	ZE-0750T
95mm	ZE-0751S	ZE-0751T




Head Dia :8 mm Hexagonal socket : 3.5mm

Core Dia : 3 mm Tap Dia : N/A

TPI Dia : 1.75 mm Drill Bit : 3.2 mm

Shaft Dia : 3 mm Guide Wire : N/A

Thread Dia: 4.5 mm



Shop No1, St#7, Main Bazar Asad Colony, Sheikhupura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com



### 3.5 mm HerBert Screw (Headless) Self Tapping

Length	Stainless Steel 316L	Titanium Grade 5
10mm	ZE-0752S	ZE-0752T
12mm	ZE-0753S	ZE-0753T
14mm	ZE-0754S	ZE-0754T
16mm	ZE-0755S	ZE-0755T
18mm	ZE-0756S	ZE-0756T
20mm	ZE-0757S	ZE-0757T
22mm	ZE-0758S	ZE-0758T
24mm	ZE-0759S	ZE-0759T
26mm	ZE-0760S	ZE-0760T
28mm	ZE-0761S	ZE-0761T
30mm	ZE-0762S	ZE-0762T



Head Dia :4.5 mm Hexagonal socket : 2mm

Core Dia : 2.7 mm Tap Dia : 3.5mm

TPI Dia : 1.25 mm Drill Bit : 2.7 mm


Shaft Dia : 2.7 mm Guide Wire : 1.2mm

Thread Dia: 3.5 mm

### ACL/PCL Screw




Length	7.0mm	8.0mm	9.0mm	10.0mm	11.0mm
20mm	ZE-0763	ZE-0768	ZE-0773	ZE-0778	ZE-0783
25mm	ZE-0764	ZE-0769	ZE-0774	ZE-0779	ZE-0784
30mm	ZE-0765	ZE-0770	ZE-0775	ZE-0780	ZE-0785
35mm	ZE-0766	ZE-0771	ZE-0776	ZE-0781	ZE-0786
40mm	ZE-0767	ZE-0772	ZE-0777	ZE-0782	ZE-0787



### DHS/DCS Screw

Length	Stainless Steel 316L	Titanium Grade 5
50mm	ZE-0763S	ZE-0763T
55mm	ZE-0764S	ZE-0764T
60mm	ZE-0765S	ZE-0765T
65mm	ZE-0766S	ZE-0766T
70mm	ZE-0767S	ZE-0767T
75mm	ZE-0768S	ZE-0768T
80mm	ZE-0769S	ZE-0769T
85mm	ZE-0770S	ZE-0770T
90mm	ZE-0771S	ZE-0771T
95mm	ZE-0772S	ZE-0772T
100mm	ZE-0773S	ZE-0773T
105mm	ZE-0774S	ZE-0774T
110mm	ZE-0775S	ZE-0775T




Head Dia :N/A Hexagonal socket : 3.5mm

Core Dia : 8 mm Tap Dia : 12.7mm

TPI Dia : 3.2 mm Drill Bit : 8 mm

Shaft Dia : 8 mm Guide Wire : 2.5mm

Thread Dia: 12.7 mm



📍 Shop No1, St#7, Main Bazar Asad Colony, Sheikhupura Road  
Gujranwala Pakistan-52250.

✉ orthozenew@yahoo.com

**Poly Axial Screw**




Length	4.0mm	4.5mm	5.0mm	5.5mm	6.0mm	6.5mm
18mm	ZE-4788PAS	ZE-4-788PAS	ZE-5788PAS	ZE-5-788PAS	ZE-6788PAS	ZE-6-788PAS
20mm	ZE-4764PAS	ZE-4-764PAS	ZE-5764PAS	ZE-5-764PAS	ZE-6764PAS	ZE-6-764PAS
22mm	ZE-4765PAS	ZE-4-765PAS	ZE-5765PAS	ZE-5-765PAS	ZE-6765PAS	ZE-6-765PAS
24mm	ZE-4766PAS	ZE-4-766PAS	ZE-5766PAS	ZE-5-766PAS	ZE-6766PAS	ZE-6-766PAS
26mm	ZE-4767PAS	ZE-4-767PAS	ZE-5767PAS	ZE-5-767PAS	ZE-6767PAS	ZE-6-767PAS
28mm	ZE-4768PAS	ZE-4-768PAS	ZE-5768PAS	ZE-5-768PAS	ZE-6768PAS	ZE-6-768PAS
30mm	ZE-4769PAS	ZE-4-769PAS	ZE-5769PAS	ZE-5-769PAS	ZE-6769PAS	ZE-6-769PAS
35mm	ZE-4770PAS	ZE-4-770PAS	ZE-5770PAS	ZE-5-770PAS	ZE-6770PAS	ZE-6-770PAS
40mm	ZE-4771PAS	ZE-4-771PAS	ZE-5771PAS	ZE-5-771PAS	ZE-6771PAS	ZE-6-771PAS
45mm	ZE-4772PAS	ZE-4-772PAS	ZE-5772PAS	ZE-5-772PAS	ZE-6772PAS	ZE-6-772PAS
50mm	ZE-4773PAS	ZE-4-773PAS	ZE-5773PAS	ZE-5-773PAS	ZE-6773PAS	ZE-6-773PAS

**Mono Axial Screw**



Length	4.0mm	4.5mm	5.0mm	5.5mm	6.0mm	6.5mm
18mm	ZE-4788MAS	ZE-4-788MAS	ZE-5788MAS	ZE-5-788MAS	ZE-6788MAS	ZE-6-788MAS
20mm	ZE-4764MAS	ZE-4-764MAS	ZE-5764MAS	ZE-5-764MAS	ZE-6764MAS	ZE-6-764MAS
25mm	ZE-4768MAS	ZE-4-768MAS	ZE-5768MAS	ZE-5-768MAS	ZE-6768MAS	ZE-6-768MAS
30mm	ZE-4769MAS	ZE-4-769MAS	ZE-5769MAS	ZE-5-769MAS	ZE-6769MAS	ZE-6-769MAS
35mm	ZE-4770MAS	ZE-4-770MAS	ZE-5770MAS	ZE-5-770MAS	ZE-6770MAS	ZE-6-770MAS
40mm	ZE-4771MAS	ZE-4-771MAS	ZE-5771MAS	ZE-5-771MAS	ZE-6771MAS	ZE-6-771MAS
45mm	ZE-4772MAS	ZE-4-772MAS	ZE-5772MAS	ZE-5-772MAS	ZE-6772MAS	ZE-6-772MAS
50mm	ZE-4773MAS	ZE-4-773MAS	ZE-5773MAS	ZE-5-773MAS	ZE-6773MAS	ZE-6-773MAS



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com

### MaxilloFacial Screw –Cross Head



Length	1.5mm SS	1.5mm TT	2.0mm SS	2.0mm TT
8mm	ZE-775MFS	ZE-786MFS	ZE-797MFS	ZE-808MFS
10mm	ZE-776MFS	ZE-787MFS	ZE-798MFS	ZE-809MFS
12mm	ZE-777MFS	ZE-788MFS	ZE-799MFS	ZE-810MFS
14mm	ZE-778MFS	ZE-789MFS	ZE-800MFS	ZE-811MFS
16mm	ZE-779MFS	ZE-790MFS	ZE-801MFS	ZE-812MFS
18mm	ZE-780MFS	ZE-791MFS	ZE-802MFS	ZE-813MFS
20mm	ZE-781MFS	ZE-792MFS	ZE-803MFS	ZE-814MFS
22mm	ZE-782MFS	ZE-793MFS	ZE-804MFS	ZE-815MFS
24mm	ZE-783MFS	ZE-794MFS	ZE-805MFS	ZE-816MFS
26mm	ZE-784MFS	ZE-795MFS	ZE-806MFS	ZE-817MFS
28mm	ZE-785MFS	ZE-796MFS	ZE-807MFS	ZE-818MFS

### 3,4mm Locking Bolt

The bolt is used for Femoral Nail



Length	Stainless Steel 316L	Titanium Grade 5
20mm	ZE-0819S	ZE-0819T
22mm	ZE-0820S	ZE-0820T
24mm	ZE-0821S	ZE-0821T
26mm	ZE-0822S	ZE-0822T
28mm	ZE-0823S	ZE-0823T
30mm	ZE-0824S	ZE-0824T
32mm	ZE-0825S	ZE-0825T
34mm	ZE-0826S	ZE-0826T
36mm	ZE-0827S	ZE-0827T
38mm	ZE-0828S	ZE-0828T
40mm	ZE-0829S	ZE-0829T
42mm	ZE-0830S	ZE-0830T
44mm	ZE-0831S	ZE-0831T
46mm	ZE-0832S	ZE-0832T
48mm	ZE-0833S	ZE-0833T
50mm	ZE-0834S	ZE-0834T
52mm	ZE-0835S	ZE-0835T
54mm	ZE-0836S	ZE-0836T
56mm	ZE-0837S	ZE-0837T
58mm	ZE-0838S	ZE-0838T
60mm	ZE-0839S	ZE-0839T

Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.


orthozenew@yahoo.com

### 4,9mm Locking Bolt

**This bolt is used for the fixation of Proximal and distal holes of pfna2, proximal femoral nail(pfn), antigrade femoral nail(AFN), Expert tibia nail, tibia nail and femoral nail**

Length	Stainless Steel 316L	Titanium Grade 5
22mm	ZE-0860S	ZE-0860T
24mm	ZE-0861S	ZE-0861T
26mm	ZE-0862S	ZE-0862T
28mm	ZE-0863S	ZE-0863T
30mm	ZE-0864S	ZE-0864T
32mm	ZE-0865S	ZE-0865T
34mm	ZE-0866S	ZE-0866T
36mm	ZE-0867S	ZE-0867T
38mm	ZE-0868S	ZE-0868T
40mm	ZE-0869S	ZE-0869T
42mm	ZE-0870S	ZE-0870T
44mm	ZE-0871S	ZE-0871T
46mm	ZE-0872S	ZE-0872T
48mm	ZE-0873S	ZE-0873T
50mm	ZE-0874S	ZE-0874T
52mm	ZE-0875S	ZE-0875T
54mm	ZE-0876S	ZE-0876T
56mm	ZE-0877S	ZE-0877T
58mm	ZE-0878S	ZE-0878T
60mm	ZE-0879S	ZE-0879T
65mm	ZE-0880S	ZE-0880T
70mm	ZE-0881S	ZE-0881T
75mm	ZE-0883S	ZE-0883T
80mm	ZE-0883S	ZE-0883T
85mm	ZE-0884S	ZE-0884T
90mm	ZE-0885S	ZE-0885T
95mm	ZE-0886S	ZE-0886T
100mm	ZE-0887S	ZE-0887T

Manufacture &  
Orthopedic



Shop No1, St#7, Main Bazar Asad Colony, Sheikhpura Road  
Gujranwala Pakistan-52250.

orthozenew@yahoo.com